

3SHAPE TUTORIAL





DIGITAL RECORDS FOR IMPLANT SUPPORTED PROSTHESIS

Complete step-by-step guide for using 3Shape intraoral camera & scanning
software for implant prosthesis

Order Setup

simplified

To begin the setup, open a new case, on the left-hand side of the screen, select the “Bridge” option. A drop-down menu will appear, select “Screw-retained.” Next, you’ll select the teeth that correspond with your implant placement.



You'll then choose “Pontic” from the drop-down menu and select the remaining teeth you want included in the prosthesis.



Implant Connection Setup

simplified

Now, proceed to Section “2. Implant Connection” This is where you’ll select the implant manufacturer, system, and connection type.

To update the implant manufacturer, click on “Manufacturer.”

A secondary drop-down menu will appear with a list of available options. Select the implant manufacturer that matches the implants being used.

The screenshot shows the 'Edit bridge' interface. On the left, a list of steps is shown: '1. Add bridge elements Screw Retained Crown,...' (checked), '2. Implant connection GM' (highlighted with a red box), '3. Material: bridge', and '4. Shade'. On the right, a secondary menu is open for '2. Implant connection', showing options: 'A. Manufacturer NEODENT' (checked), 'B. System GM' (checked), and 'C. Implant connection'. At the bottom, there are 'Save' and 'Cancel' buttons. The top right shows 'Tooth 3: screw-retained... GM' with a tooth icon.

The screenshot shows the 'Edit bridge' interface. On the left, a list of steps is shown: '1. Add bridge elements Screw Retained Crown,...' (checked), '2. Implant connection' (highlighted with a red box), '3. Material: bridge', and '4. Shade'. On the right, a secondary menu is open for '2. Implant connection', showing options: 'A. Manufacturer Straumann' (checked and highlighted with a red box), 'NEODENT', 'Nobel Biocare', and 'Straumann' (highlighted with a red box). At the bottom, there are 'Save' and 'Cancel' buttons. The top right shows 'Tooth 3: screw-retained...' with a tooth icon.

After selecting the manufacturer, you’ll be prompted to choose the specific implant system and connection from the options provided.

This process will need to be repeated for each implant site to ensure accurate selection.

For quicker future selections, you can set up “favorites” within your system.

Material Setup

simplified

Now, proceed to Section “3. Material” In this section, you’ll select the material type and, if applicable, the desired shade.

Otherwise, select “I don’t want to set a shade,” then save your order data to complete this step.

This screenshot shows the 'Edit bridge' screen with a list of steps on the left: 1. Add bridge elements (Screw Retained Crown, PMMA), 2. Implant connection (Bone Level, RC), 3. Material: bridge (PMMA), and 4. Shade. The third step is highlighted with a red box. On the right, a list of material options is shown: Emax, Gold, PMMA (highlighted with a red box), Zirconia - monolithic, and Zirconia - veneered. A 'Save' button is at the bottom left.

This screenshot shows the 'Tooth 3..14' screen. On the left, a list of steps is shown: 1. Add bridge elements, 2. Implant connection, 3. Material: bridge (PMMA), and 4. Shade. The fourth step is highlighted with a red box. On the right, a list of shade options is shown: A1, A2, A3, A3.5, A4, B1, B2, and B3. A red box highlights the option 'I don't want to set a shade' at the top. A 'Save' button is at the bottom left.

This screenshot shows the 'Edit bridge' screen with the same list of steps as the previous screenshots. The 'Save' button at the bottom left is highlighted with a red box.

This screenshot shows the 'Patient' screen. Under the 'Case overview' section, there are two items: 'Bridge 3..14' and 'Screw Retained Crown 3, 5, 7, 10, 12, 14'. Under the 'Additional scans' section, the 'Pre-preparation' option is checked and highlighted with a red box. Other options include 'Emergence profile' (checked) and 'Enable intraoral camera' (unchecked).

On the right-hand side of the screen, make sure to check the “Pre-Preparation” option under Additional scans. This will allow you to scan the teeth, gingiva, and scan bodies all within a single case order.



Scanner Settings

simplified

When scanning the teeth intraorally, it's best to keep the AI feature turned ON to ensure accurate capture of the tooth surfaces.


However, when scanning the interface between the bridge and the patient's gingiva, you'll want to turn the AI OFF. This allows the scanner to better capture the soft tissue detail without AI removing valuable anatomical data.

To disable the AI feature, simply tap the button at the bottom of the screen while the scanner is active.



Alternatively, you can turn off the AI setting by navigating to the "Tools" menu on the left-hand side of the screen.

Once inside the Tools menu, simply uncheck the "AI Scan" option to disable it.






Scanning Prostheses

simplified

You will now begin scanning. Always start by scanning the temporary teeth. When doing so, ensure you capture enough of the surrounding gingiva for the scans of the teeth and gingiva to properly mesh.

Be sure to capture the palate, and when scanning the mandibular arch try to capture the ridge behind the posterior teeth. This tends to provide the best results when the software meshes the lower scans.


Do your best to clearly capture the junction between the prosthesis and the gingiva.



Pro tip: Capture the bite before proceeding with the gingiva and scan body scans. The software will prompt you to take both the left and right bite.

Gingival Scan


simplified



Before starting the gingival scan, the software will prompt you to mark the implant locations.

Once marked, the selected areas will appear highlighted in blue, each implant site will be confirmed with a green checkmark at the bottom of the screen.

Press "Done" to proceed with the gingival scan.




You'll notice that the previously blue-marked areas will be cleared, allowing for a more accurate capture of the gingival tissue.

Once you've scanned the entire arch, the teeth will disappear from view, leaving you with a scan that displays only the gingiva for precise visualization and design reference.



Scanbody Scan

simplified





You'll follow the same steps as you did prior to the gingival scan: begin by marking each implant site, then press "Done" to proceed.

Next, carefully place the appropriate scanbodies onto the implants, following the manufacturer's guidelines. Ensure each scanbody is fully seated and correctly aligned.

Once the scanbodies are in place, you'll begin scanning. Make sure to capture all relevant surfaces—including the top, sides, and surrounding gingiva for precise alignment and accurate implant positioning.

Before finalizing the scan, double-check that all scanbodies are clearly visible and have been fully captured.




Verified Model Scans


Once these scans are complete, you'll need to open a new case to scan the arches on your verified models. Ensure the models include indicators, as these are essential for accurately aligning any future model scans performed at the lab with the intraoral scans you've taken.

When scanning the teeth on the models, be sure to enable "Model Scanning" under the Tools menu. This setting helps ensure optimal accuracy and detail during the scanning process.

Do your best to capture the joint between the analogs and the arches as clearly as possible.



simplified



“At Simplified, our mission is to transform complex dental processes into streamlined workflows, delivering exquisite prosthetics crafted from the highest quality materials. We are committed to setting realistic expectations and consistently exceeding them in thoughtful, innovative ways—ensuring both clinicians and patients experience exceptional results at every step.”

info@simplified123.com

8043 Wornall Rd, Ste 102
Kansas City, MO 64114

www.simplified123.com