

# SHINING3D TUTORIAL

## DIGITAL RECORDS FOR IMPLANT SUPPORTED PROSTHESIS

Complete step-by-step guide for using Shining3D intraoral scanning software  
for implant prosthesis

---



# New Order

simplified

To enter a new order, select the New Order button located at the top left of your screen.

If taking scans on an existing patient within the software, enter the Patient's name in the "Patient Name" section & select the corresponding Patient.

If adding a new patient into the software select the add patient button – see below:



# Adding a New Patient

simplified

The screenshot displays the SHINING3D software interface. At the top, a progress bar shows four stages: 'New Order' (highlighted in blue with 'In Progress' status), 'Scan', 'Pre-design', and 'Send'. The 'New Order' section on the left contains fields for 'ID' (095), 'Create Time' (2/18/25 10:27 AM), 'Order Number' (095), 'Type' (First Visit), 'Patient Name\*' (empty), 'Doctor' (Doctor001), 'Operator' (001 Technician001), 'Lab' (General), and 'Dentistry Type' (Implant selected). A 'Notes' section is at the bottom left. The 'New Patient' pop-up form is centered, featuring fields for 'Name\*', 'Gender\*' (Male), 'Age\*', 'Date of Birth', 'Contact', 'Pregnant' (No), 'Institution' (Simplified LLC), 'Medical History' (Diabetes, Hypertension), 'Allergies', and 'Oral Habits' (Finger Biting, Lip Biting, Bruxism, Oral Breathing, Smoking, Sugar). 'Confirm' and 'Cancel' buttons are at the bottom of the pop-up. The background interface includes buttons for 'Full Crown', 'Bridge', 'Upper Jaw', 'Lower Jaw', 'Full Jaw', and a 'Scan' button at the bottom right.

After selecting “add patient” a pop-up will appear.

Enter the Patient information—anything that has an asterisk next to it is required.


Once all the information is entered click the confirm button at the bottom of the pop-up window.



# Kit Selection Process

simplified

1. Towards the bottom & center of the screen you'll see the following:




2. The serial number in orange indicates the scanning kit that you will be utilizing, you can find the serial number on the kit here:




\*\*these scan flags need to be put back in the correct kit, the serial number is located on the underside of each scan body:




3. To change the kit click the button located directly next to the serial number within the software.



4. This will prompt another pop-up window with a drop-down menu to select the correct serial number that corresponds to whichever scan body kit you are using.



5. After selecting the correct serial number you will then select Upper Arch or Lower Arch or both.



# Arch Selection Setup



simplified

Next, you'll select the abutment locations—this step must be accurate. If the locations are incorrect, the software will flag the issue during scanning and prevent you from proceeding. Additionally, correcting the abutment locations mid-scan will delete all scans taken up to that point.

Be sure to also indicate the correct number of implants before moving forward.

Once you've finished selecting the abutment locations, be sure to select the corresponding arches under the "Scan a Pre-op Model" section on the right-hand side of the screen.

After saving the information, proceed by clicking the Scan button. This will prompt you to select the implant system—click Yes to continue.



A screenshot of a software settings panel titled 'Scan A pre-op Model'. The panel is outlined with a red border. It contains the following sections: 'Implant systems' with a 'Select Implant >' link; a note '\*Please select the tooth position first'; checkboxes for 'Upper Jaw' and 'Lower Jaw', both of which are checked; an 'Implant Base' dropdown menu set to 'Custom Abutment'; a 'Material' dropdown menu set to 'Composite Material'; and a 'Tooth Shade' dropdown menu set to 'None'. At the bottom of the panel are three buttons: 'Explore', 'Save', and 'Scan'.



# Implant System Selection

simplified

A new pop-up window will appear displaying the implant systems the office frequently uses or has used in the past.



When scanning for full arch cases always choose the 4mm Gold from ImplantDirect. This is a DESS scan body that is universal and the majority of labs will request. Select the "All Applications" button and then hit the confirm button.


# Scanning Prostheses

simplified

You will now begin scanning. Always start by scanning the temporary teeth. When doing so, ensure you capture enough surrounding gingiva so that the scans of the teeth and gingiva can properly mesh.

For the maxillary arch, be sure to capture the palate, and for the mandibular arch, try to capture the ridge behind the posterior teeth—this tends to provide better results when the software meshes the lower scans.

Also, do your best to clearly capture the junction between the prosthesis and the gingiva.



Pro tip: Capture the bite before proceeding with the gingiva and scan body scans. The software will prompt you to take both the left and right bite.

As long as you receive green checkmarks next to each scan, you are ready to move forward with the tissue scans.




# Scanning Gingiva

simplified

Once you've verified the abutments are torqued to the manufacturers recommendation and completed the jig verification, you can proceed with the gingiva scans. The use of an Optragate is recommended to help obtain clear and consistent scans.

Be sure to capture the palate on the maxilla, as well as the buccal/facial angles on both the maxilla and mandible for optimal accuracy.



If there is adequate data, the initial scan with the temporary teeth in place and the gingiva scan will automatically mesh. If the scans do not mesh, it typically means that there wasn't enough information captured—either in the teeth scan or the gingiva scan.



When the scans mesh correctly, you'll see two green checkmarks at the bottom of the screen, and the result will look like the example shown on the left.



# Scanning Flags

simplified

You'll now place the scan flags. The screws within the flags are Straumann and should be torqued to 10 Ncm. Ensure that the flags do not overlap or interfere with each other's ability to seat fully. Recommendation; use an explorer to trace around the margin of each flag and abutment to confirm they've been fully seated. When selecting flags, use the longer scan flags in the posterior and the shorter ones in the anterior for optimal visibility and accuracy.




After the software recognizes all scan flags, it will prompt you to scan the junction where each flag meets the abutment—on both the right and left sides. When the scans mesh correctly, you'll receive four green checkmarks at the bottom of the screen.

# Completing Your Scans

simplified


When moving on to the second arch or finishing your scans, the software will again prompt you to select the implant system. Always choose the 4mm Gold from ImplantDirect.

At this point, your scans should be complete and ready to send or export.





# simplified



“At Simplified, our mission is to transform complex dental processes into streamlined workflows, delivering exquisite prosthetics crafted from the highest quality materials. We are committed to setting realistic expectations and consistently exceeding them in thoughtful, innovative ways—ensuring both clinicians and patients experience exceptional results at every step.”

[info@simplified123.com](mailto:info@simplified123.com)

8043 Wornall Rd, Ste 102  
Kansas City, MO 64114

[www.simplified123.com](http://www.simplified123.com)